

DAWN INR CAPTURE

A quick and easy, web-based, mobile / tablet friendly way for INR self-testers, domiciliary services and remote phlebotomy services to capture INR results and supporting information and send it to physicians, specialist anticoagulation therapy management systems such as DAWN AC or electronic medical records.

In-built checks and verification mean that the INR result can be imported straight into the patient's DAWN AC record.

DAWN INR Capture is compatible for use with PC, Laptop, Tablet and Phone.

Example workflow:

STEP 1

The patient is issued with a PIN number generated by and recorded in DAWN AC and issued via SMS, email or letter.

STEP 2


When the patient self-tests, they complete the DAWN INR Capture form (using the PIN and DOB to identify themselves - see image on the right) and clicks 'post answers'.

STEP 3

DAWN INR Capture sends the results to the DAWN AC patient record where dose and next test date are calculated.

STEP 4

In range patients are sent dosing instructions via SMS, email or letter. Out of range patients are contacted by the healthcare professional. (*This element of the workflow depends on local processes and procedures*).



Record your INR.

Please answer all questions

PIN

DOB

Date

INR

Have you missed any warfarin doses? Yes No

Have you taken too many warfarin tablets? No

Have you experienced any bleeding or bruising? No

Have you had any change to your medications? No

Have you started taking any non-prescription or herbal medications? No

Have you been admitted to hospital or had any surgical procedure? No

Please enter details for any Yes answers

Missed dose on Wed

Post Answers

