DawnMailer Troubleshooting

DawnMailer Troubleshooting	1
Setting Up DawnMailer: checklist	2
Run-time Error 380 – invalid property value	5
Error: -2146697210, the system cannot locate the object specified	7
Error: no or invalid MessagesInterfaceKey given	8
Problem: error message still appears at the bottom of the screen	9
Question: how can I tell if DawnMailer is configured correctly?	10
Question: do I need to press the Submit Printers button?	10
Question: why do some of my printers now appear as Not In Use?	11
Problem: message failed to print	12
Problem: emails not received	18
Undeliverable Emails: general troubleshooting	19
Undeliverable Email: the server rejected the sender address	20
Undeliverable Email: the transport failed to connect to the server	21
Diagram: Print Station / Dawn Mailer Printer Mappings	23

Setting Up DawnMailer: checklist

What do I need to check before I set up DawnMailer?

Check #1

Choose Start/Printers and Faxes from the Windows Start menu (on the DawnMailer machine) to display the printers.



- 1. Ensure all local and/or network printers you need to print to appear in the list.
- 2. Print a test page (or ideally a web page using Internet Explorer) to each printer to ensure it prints correctly.
- 3. Ensure any printing problems you encounter at this stage are resolved before you attempt to configure DawnMailer.

Check #2

If DawnMailer is to run as a service, it is preferable to have a specified user login it can run under, rather than running under the Local System account. DawnMailer uses Internet Explorer to print messages. Internet Explorer creates a Temporary Internet File for each page it opens. While this boosts efficiency for normal web browsing, it does not provide any benefit for Dawn CF / DawnMailer but it can use a very large amount of disk space.

If DawnMailer runs under the local system account, the Temporary Internet Files folder is created under the Documents and Settings folder for the Default User. This is usually on the c:\ drive (often reserved for the operating system on Windows Servers) and can create performance issues for some remote access technologies or virus checkers that routinely scan this folder.

If DawnMailer runs under a specified user login, you can log on to the server as this user and set their Internet Explorer Preferences to limit the size of the Temporary Internet Files folder. If necessary, you can also change the location of the folder to a different drive.

Ensure you have a local or domain user account, which you can run use to run DawnMailer. It is **imperative** that the **password** for this account is configured to **never expire**. If this is not possible, you will have to run DawnMailer under the local system account.

Once you have the account's login details:

- 1. Login to the server as this user and open Internet Explorer.
- 2. Choose Tools/Internet Options from the menu.
- Click the Settings button under *Temporary* Internet Files or Browsing History (depending on the browser version) to display the Temporary Internet Files Settings.
- 4. Set the *Disk Space to Use* setting to the lowest allowable value (this may vary slightly depending on your operating system and browser version; 1 Mb or 8 Mb are common).
- 5. If necessary use the Move folder button to change the location of the Temporary Internet Files folder for this user.

Te	mporary Internet Files and History Settings 🛛 🛛 🔀				
ſ	Temporary Internet Files Internet Explorer stores copies of webpages, images, and media for faster viewing later.				
	Check for newer versions of stored pages: O Every time I visit the webpage				
	O Every time I start Internet Explorer				
	<u> </u>				
	Disk space to use (8 - 1024MB): 8 📚				
	Current location: C:\Documents and Settings\George.4SLIVE\Local Settings\Temporary Internet Files\				
	Move folder View objects View files				
ſ	History				
	Specify how many days Internet Explorer should save the list of websites you have visited.				
	Days to keep pages in history: 5 💲				
	OK Cancel				

Check #3

Check your login (or the login under which DawnMailer will run if different) has Write permissions for the folder in which DawnMailer resides and all of its subfolders.

To check this:

- 1. Open the folder in Windows Explorer, right click on any bit of empty white space and choose New/Text Document from the pop up context menu. A new Text document is created. Double click the new text document, add a line of text and choose File/Save.
- 2. If you can create and update the file in this way you have sufficient permissions.
- 3. If you receive an "Access Denied" error message you do not have sufficient permissions.
- 4. DO NOT attempt to run DawnMailer until you have resolved the permissions shortfall (you will need to ask IT to do this as if you do not have Write Permissions, you almost certainly will not have permission to alter Folder Permissions).

Run-time Error 380 – invalid property value.

Run-time error 380 – invalid property value. I can't open DawnMailer. When I try to do so, I get this error message. When I click OK on the message, DawnMailer fails to open.

Solution

See *Setting Up DawnMailer: Checklist – Check #1* above: check you have at least one printer installed on the machine on which you are trying to run DawnMailer.

If IT have installed the printers but you cannot see them if you choose *Start/Settings/Printers & Faxes*, report this to IT. They may have installed the printers in such a way that they are only visible to the log in they were installed under. They should be available to everyone.

If you had printers installed on the machine and DawnMailer has been working but you are now experiencing this error, choose *Start/Settings/Printers & Faxes* to make sure the printers are still visible. If they are not check the Print Spooler service is running.

To check the Print Spooler service is running:

- 1. Right click on My Computer and choose Manage from the pop up context menu to display the Computer Management Console.
- 2. Choose Applications & Services/Services to show the list of services on the machine.
- 3. Locate Print Spooler in the list and ensure its status is *Started*.
- 4. If not right click on the Print Spooler service and choose Properties from the context menu.



5. Ensure the Startup Type option is set to *Automatic* and choose the Start button.

Print Spooler Pr	operties (Local Computer)	?×			
General Log On	Recovery Dependencies				
Service name:	Spooler				
Display <u>n</u> ame:	Print Spooler				
<u>D</u> escription:	Loads files to memory for later printing.	~			
Pat <u>h</u> to executab	le:				
C:\WINDOWS\s	ystem32\spoolsv.exe				
Startup type: Automatic					
Service status:	Stopped	_			
<u>S</u> tart	Stop Pause Resume				
You can specify from here.	the start parameters that apply when you start the serv	/ice			
Start parameters:					
	OK Cancel	spply			

6. Choose OK to save the changes.

Error: -2146697210, the system cannot locate the object specified

Error: -2146697210, The system cannot locate the object specified

...I get the above error message when I start DawnMailer

Solution

Check you have entered the correct server address in the DawnMailer configuration screen. The server address is the web address (URL) or your Dawn application WITHOUT index.html on the end, e.g.

http://DawnServer/DawnAC

...not

http://DawnServer/DawnAC/Index.html

📾 DawnMailer 1.0 (27)			
Server address	http://localhost/Dawr	AC/	
Use CDOSYS (recommended)	🔽 Halt on	errors 🕅	
CDOSYS SMTP server			
Default printer	\\4SDC1\hp LaserJet	1320 PCL 6 💌 📃	Submit
Interval (sec.)	5 💌]	
Interface key	9AF62AC2E0D2B72A		Test Mail Test Email
		_	Test Fax
🗆 # messages sent			
Total	0 Printed		0
	Email		0
Time of last message	Fax	Vie	0 wlog l
			milog
			Statt)

Error: no or invalid MessagesInterfaceKey given

Error: No or invalid MessagesInterfaceKey given

...I get the above error message when I start DawnMailer

Solution

Copy and Paste correct Identification Key from your Print Station in Dawn AC into the InterfaceKey in the DawnMailer configuration screen.

📾 DawnMailer 1.0 (27)				_ 🗆 🗙
Server address	http://loc	alhost/DawnAC/		
Use CDOSYS (recommended)		Halt on errors 🔲		
CDOSYS SMTP server				
Default printer	\\4SDC1	hp LaserJet 1320 PCL 6	-	Submit
Interval (sec.)	5	•		
Interface key	SAF62AC	2E0D2B72A		Test Mail Test Email
				Test Fax
# messages sent				_
lotal	U	Printed		U
		Email		U
Time of last message		rax		View log
				Ticknog
				Start

Problem: error message still appears at the bottom of the screen

I have corrected the problem but the error message still appears at the bottom of the screen.

Solution

DawnMailer does not clear the last error message automatically. If a new error occurs, it is displayed. If no further error occurs, the last error remains on screen until you closed and reopen DawnMailer. Close and Reopen DawnMailer to clear the last error message.

📾 DawnMailer 1.0 (27)				
Server address	http://local/	nost/DawnAC/		
Use CDOSYS (recommended)		Halt on errors 🔲		
CDOSYS SMTP server				
Default printer	\\4SDC1\h	p LaserJet 1320 PCL 6	•	Submit
Interval (sec.)	5	-		
Interface key	9AF62AC2E	:0D2B72A		Test Mail
				Test Fax
# messages sent				
Total	0 F	Printed		0
	E	Email		0
	F	Fax		0
Time of last message -				View log
	. ·]		
No or invalid MessagesInterfac	eKey given			Stop

Question: how can I tell if DawnMailer is configured correctly?

How can I tell if DawnMailer is configured correctly to connect to Dawn AC?

Answer

Open your Print Station in Dawn and check the list of printers. When you start DawnMailer, it submits the list of printers it can see to Dawn, which stores them in the database under the Print Station. If DawnMailer is correctly configured, the printers now appear under the Print Station.

Printers			
			~
Name	DeviceName	PrinterPaperType In	U
資 \\4SDC1\hp LaserJet 1320 PCL 6	\\4SDC1\hp LaserJet 1320 PCL 6	×	
\\4sdc1\Colour	\\4sdc1\Colour	X	1
\\4SDC1\HP CLJ2820 PCL 6 - Black_White	\\4SDC1\HP CLJ2820 PCL 6 - Black_White	X	
Microsoft Office Document Image Writer	Microsoft Office Document Image Writer	X	1
FinePrint	FinePrint	X	<u>, </u>
🗑 Fax	Fax	X	1
🖎 Adobe PDF	Adobe PDF		
<		>	
<< < > >>	New Delete Save	To list Pri	int

Please note: if your Print Station screen is already open, you must close it and reopen it to refresh the contents and display the printers.

Question: do I need to press the Submit Printers button?

Do I need to press the Submit Printers button to send the list of printers to Dawn?

Answer

No. When you start DawnMailer, it automatically submits the printers to Dawn. Use the Submit Printers button to resubmit the printers (if you have made any changes), while DawnMailer is running or while it is in Stopped mode.

Question: why do some of my printers now appear as Not In Use?

Since the last time I started DawnMailer, some of my printers now appear as Not In Use in the Print Station List.

Answer

DawnMailer can no longer see all the same printers it could see before. You may have:

- Uninstalled certain Windows printers from the DawnMailer machine
- Moved DawnMailer to another machine
- Changed the login DawnMailer is running under to one that cannot see the same printers as before

Check the available Windows printers on the DawnMailer machine (see Question #1 above). These match the printers marked as In Use in the Print Station. DawnMailer automatically marks any Print Station printers that it can no longer see as not In Use

Printers				
		D.:	()	^
Name Name \\4SDC1\hp LaserJet 1320 PCL 6	\\4SDC1\hp LaserJet 1320 PCL 6	PrinterPape	riype Inuse	
\\4sdc1\Colour	\\4sdc1\Colour			_
\\4SDC1\HP CLJ2820 PCL 6 - Black_White	\\4SDC1\HP CLJ2820 PCL 6 - Black_White		X	=
🚔 Microsoft Office Document Image Writer	Microsoft Office Document Image Writer			
🚔 FinePrint	FinePrint		X	
😭 Fax	Fax		X	
🖏 Adobe PDF	Adobe PDF		X	⊻
<u><</u>			>	

DawnMailer did not print my message.

Initial Troubleshooting

- 1. Check the Message Recipient Status for the Message by performing the following steps:
- 2. Click the message in the List of messages on the Patient's Letters Tab to display the Message.
- 3. Click the Recipients Tab to display the recipient or recipients and check the status column.

1	Message	Recipients Workflow			
	Deliv	veryDate RecipientHCprofessional	RecipientPatient Sent by	Status	LastStateChange
	9	Boots The Chemist - Pinter John	Email	Sent	15/07/2008 16:18

Status	Meaning	What to check
Sent	The message has been sent into the Dawn Message queue but has not yet been processed by DawnMailer. Messages with a status of "Sent" are often referred to as "Pending" messages	 Close the Recipients form. Wait a few seconds and reopen it to see if the status has changed. If Dawn is busy processing many other messages, it may be slower than normal in printing yours. If the status does not change, DawnMailer may have stop running. Restart DawnMailer. See Solution #1 below to check DawnMailer is configured to run unattended
Delivered	DawnMailer has delivered the message to the printer.	Choose Printers & Faxes from the Windows Start Menu on the DawnMailer machine. Double Click the appropriate printer to view the contents of its printer queue. If the Message is in the queue, check that the document at the top of the queue has not stopped printing with an error (out of paper, paper jam etc). If so, resolve the printer error. If the printer queue is empty or does not contain your message, see Solutions #2 & #3 below to check the Printer mappings are correct and HTMLPrinter WaitTime setting is adequate.
Undeliver- able	DawnMailer encountered an error trying to send the message to the printer	See Solution #4 below on dealing with Undeliverable Messages

Solution #1 - Sent (i.e. pending) Messages

DawnMailer connects to the Dawn application every 5 seconds to check if there are any pending messages. However, if Dawn is busy (perhaps processing a bulk message for a long list of recipients) it may not be able to respond to DawnMailer immediately. In this case, DawnMailer thinks an error has occurred and may switch itself into Stopped mode. If DawnMailer is running unattended on a server or on a PC that is not readily accessible by the users, this can be inconvenient especially where trying again in another 5 seconds would have resolved the problem.

In DawnMailer version 1.0.0.21 or later, you can configure DawnMailer not to Halt On Errors. If so configured, DawnMailer does not switch itself into stopped mode when it encounters an error. Instead it writes the error to a date and time stamped file in an Errors sub folder then tries again in another 5 seconds.

If DawnMailer stops processing messages and simply restarting DawnMailer fixes the problem, consider configuring DawnMailer not to Halt On Errors, by performing the following steps:

1. In order to make this change you need to run DawnMailer as an application rather than a service. If it is already running as an application, go straight to step 2 below. However, if DawnMailer is running as a service under FireDaemon, first open FireDaemon and stop the service by highlighting DawnMailer and pressing the Stop video button.

🛛 FireDaemon OEM Service Manager v1.9 GA	X
<u>File S</u> ervice <u>H</u> elp	
Service Status Process Startup Type User	
🗑 FireDaemon Service: Dawn Mailer Running Running Automatic Local System	
[] [>

- 2. Double click the DawnMailer.exe (usually located in the Polar Studio\Webroot folder) to display the DawnMailer configuration form.
- 3. If the fields are greyed out, press the Stop button.

4. Uncheck the Halt On Errors checkbox.

📾 DawnMailer 1.0 (27)			
Server address	http://loc	alhost/DawnAC/	
Use CDOSYS (recommended)	$\overline{\mathbf{v}}$	Halt on errors 🔲	
CDOSYS SMTP server			
Default printer	\\4SDC1	hp LaserJet 1320 PCL 6	➡ Submit
Interval (sec.)	5	-	
Interface key	9AF62AC	2E0D2B72A	Test Mail
			Test Fax
# messages sent			
Total	29	Printed	1
		Email	2
Time of last message 17/0	17/2008 13	:16:05	View log
			(Statt)

If you do not have the Halt on Errors checkbox, you are using an earlier version of DawnMailer and you will need to upgrade to the latest version of Dawn before you can take advantage of this feature.

(If you are running version 7.6, you may be able download version 1.0.0.21 of DawnMailer from the 4S Distrib FTP Site – contact 4S for instructions. DawnMailer version 1.0.0.21 works with version 7.6 but did not originally ship with it).

5. Press the Start button so DawnMailer saves its new configuration.

If DawnMailer normally runs as a Service, you must stop DawnMailer running as an application (as it is now) and restart it as a service, by performing the following steps:

- 6. Press the Stop button to Stop DawnMailer running as an application (if it has started processing messages, it may not respond to the Stop button until it has finished). Then close the DawnMailer configuration window.
- 7. Restart DawnMailer under FireDaemon by pressing the Start video button.

Solution #2 – Delivered messages that do not appear to have printed

Check whether the messages have printed on a different printer to the one you were expecting.

If so, you need to set up your Paper Types and Locations in Dawn (see the *Setting Up Printing – User Guide* section of the Dawn AC Help File, shipped with v7.7 +, or contact 4S for a How To... instruction sheet).

If you want all messages to print to a single printer, ensure this printer has the lowest Order Number in the Dawn Print Station.

Name \\\4SDC1\\hp LaserJet 1320 PCL 6 \\\4sdc1\Colour \\4SDC1\\HP CLJ2820 PCL 6 - Black_White Microsoft Office Document Image Writer FinePrint	DeviceName \\4SDC1\hp LaserJet 1320 PCL 6 \\4sdc1\Colour \\4SDC1\HP CLJ2820 PCL 6 - Black_White Microsoft Office Document Image Writer FinePrint	PrinterPaperT	
Name \\4SDC1\hp La	aserJet 1320 PCL 6		
DeviceName \\\4SDC1\hpLa	aserJet 1320 PCL 6		
Notes			
PrinterPaperType (None selecte	ed)		
OrderNr 1			
InUse L PrinterLocation Interest of the second seco	tems to display		

If your messages are marked as Delivered and you are quite certain they have not printed on any of the printers listed in the Print Station, see solution #3 below.

Solution #3 – Delivered Messages that have not actually printed

When DawnMailer receives a message to print, it writes the contents to a file called Temp.HTML. It then uses another program, HTMLPrinter, to perform the printing. HTMLPrinter opens the Temp.html file in an instance of Internet Explorer, then calls the Explorer File/Print command to print the message to the designated printer.

This method ensures Dawn is always using an industry standard product to print the HTML content of the message. However, Internet Explorer does not have a means of telling HTMLPrinter when it has finished printing. As such, HTMLPrinter waits for a specified number of seconds before it closes Internet Explorer and moves on to the next message.

The WaitTime setting within the HTMLPrinter.ini file specifies the number of seconds to wait. While the default setting (currently 5 seconds) should work with the majority of printers, printers and networks vary so some may take longer than others to spool print jobs. If HTMLPrinter closes Internet Explorer before the print job has finished, the print job is lost and the message does not print (even though DawnMailer thinks it has).

In this case:

- 1. Locate the HTMLPrinter.ini file in the DawnMailer folder.
- 2. Stop DawnMailer
- 3. Open this file using Notepad or any text editor.
- 4. Increase the number of seconds specified for the WaitTime, e.g.

WaitTime = 7

5. Restart DawnMailer and check messages are now printing.

Tip: if your new time is still too small, the same problem will occur. Start by specifying a long WaitTime such as 15 seconds, to verify this does actually fix the problem. Then experiment by shortening the time until the problem reoccurs. The optimum value is probably two or three of seconds longer than the longest setting, which causes the problem.

If you make the setting too long, you introduce an unnecessary delay between the printing of each message.

Solution #4 – Undeliverable Messages

Because most printers spool print jobs to a queue and then print them in the background, some common printing errors (such as paper jams or out of paper errors) occur after DawnMailer has processed the message. DawnMailer does not know about these errors so these messages are marked as Delivered (see Delivered Messages in the Initial Troubleshooting section above).

However, if DawnMailer encounters an error while its is processing the message, it marks the message as "Undeliverable" and writes the error details to the Delivery Notes field for the Message Recipient.

Date sent Recipient Sent by Status	17/07/2008 13:15 Patient Mail Undeliverable Re-send Discard
Message info Delivery	History
Error while printi printer called "\\	ng document: -2147352567, There is no SONYVAIO2GHZ\Smart Label Printer 240".

To find out why a message is undeliverable:

- 1. Click the Message in the Patient's Letters Tab to display the Message Details
- 2. Click the Recipients tab to display the Message Recipient or Recipients.
- 3. Click the Recipient whose status is marked as Undeliverable and select the Delivery tab to see the Delivery Notes.

The resolution will depend on the reason for failure. The above example shows DawnMailer could not find the Smart Label Printer it was configured to use. The printer had perhaps been removed since DawnMailer was last started. As a consequence the printer list in the Print Station is out of date.

In this case, you should restart DawnMailer (or press the Submit printers button) to update the Dawn Print Station with the current list of printers. This marks any printers that are not longer available as not In Use. You would either then designate another printer as the label printer or reinstall the Smart Label printer and restart (or resubmit printers in) DawnMailer to update the printer list again.

Once you have resolved the problem either:

Press Re-Send to re-send the message to the printer or Discard to remove the message from the undeliverable list.

Emails are marked as Delivered but never arrive.

Solution

Check the email address for the patient or healthcare professional is correct. If it is not correct, the email is sent out but never arrives. An undeliverable email notification may be sent back to the sender email address but this is not logged in Dawn.

If the email address is correct check your DawnMailer settings as follows:

1. If DawnMailer is running as a service under FireDaemon, open FireDaemon and stop the service by highlighting DawnMailer and pressing the Stop video button.

G FireDaemon OEM Service	e Manager v1.9 GA			
<u>F</u> ile <u>S</u> ervice <u>H</u> elp				
"≞ ≌ × × ▸	Stop the selected service	2 🔽 🗣		
Service	Status Process	Startup Type	User	
🔞 FireDaemon Service: DawnMai	iler Running Running	Automatic	LocalSystem	
<				>

- 2. Double click the DawnMailer.exe (usually located in the Polar Studio\Webroot folder) to display the DawnMailer configuration form.
- 3. If the fields are greyed out, press the Stop button.
- 4. Ensure the Use CDOSYS check box is checked and the CDOSYS SMTP Server contains the address of your SMTP Server.

If you do not enter these details DawnMailer attempts to write the email to a pick up folder on your web server for a virtual SMTP server to handle. If the virtual SMTP service is not enabled or correctly configured, the email is not delivered, but DawnMailer does not know so cannot mark the message as Undeliverable

📾 DawnMailer 1.0 (27)				
Server address	http://loc	alhost/DawnAC/		
Use CDOSYS (recommended)		Halt on errors 🔲		
CDOSYS SMTP server	smtp.MyS	MTPServer.com		
Default printer	\\4SDC1	\hp LaserJet 1320 PCL 6	•	Submit
Interval (sec.)	5	-		
Interface key	9AF62AC	2E0D2B72A		Test Mail
				Test Fax
🗆 # messages sent				
Total	60	Printed		4
		Email		3
		Fax		0
Time of last message 18/0)7/2008 15	i:20:43		View log
				Start

My emails are all being marked as Undeliverable

Solution

Find out why a message is undeliverable by checking the delivery notes as follows:

- 1. Click the Message in the Patient's Letters Tab to display the Message Details
- 2. Click the Recipients tab to display the Message Recipient or Recipients.
- 3. Click the Recipient whose status is marked as *Undeliverable* and select the Delivery tab to see the Delivery Notes.



The solution depends on the error message. Some common problems are listed below:

Undeliverable Email: the server rejected the sender address

Undeliverable email -Delivery Notes say:

An error occurred: -2147220978, The server rejected the sender address. The server response was: 550 <noreply@medical.org>... Not permitted

Solution

The sender email address you have entered in Dawn AC is not recognised or permitted by your SMTP Server. You may still be using the default entry - <u>noreply@medical.org</u>.

Change the Sender email address in Dawn to a valid email address that your SMTP server recognises and permits. You can set the default Sender email address for all organisations by editing the *Message_DefaultEmailFromAddress* System Setting. You can override this for a particular organisation by entering an organisation specific email address against the Organisation within the Print Station

Supported organisations	
Organisation ReturnEmailAd	dress
HOLYWELL SURGERY	
Printers	
Name	DeviceName \\4SDC1\hp Laser
\\4sdc1\Colour \\4sdc1\Colour \\4sbc1\BB CI 12820 PCL 6 - Black White	\\4sdc1\Colour
Microsoft Office Document Image Writer	Microsoft Office [
<< > >>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>	New

Undeliverable Email: the transport failed to connect to the server

Undeliverable email -Delivery Notes say: An error occurred: -2147220973, The transport failed to connect to the server.

Solution

This is a general message that indicates that DawnMailer failed to connect to the SMTP server you specified. This may be because:

- You have entered an incorrect SMTP server address
- A firewall on your DawnMailer machine is blocking access
- The SMTP server does not allow relaying of emails from authorised machines

... or some other reason.

Check you have entered the SMTP server address correctly.

Check with your IT department that no firewall software is blocking the SMTP port (port 25) on the DawnMailer machine.

Check with your SMTP server administrator that your DawnMailer machine is permitted to relay emails via the SMTP server. If not, ask if the DawnMailer machine can be added to a permitted list.

Alternatively, you can use Telnet commands to try and troubleshoot the problem further. The following procedure describes how to attempt to connect to the SMTP server and send a message via a Telnet session. This will not solve the problem but it may provide more precise error messages, which may help you pinpoint the source of the problem.

Please take care to type the commands exactly as shown including colons. You cannot use the backspace key to correct them if you make a mistake (backspace appears to work but the command will fail). In this example, we'll assume your SMTP server address is smtp.myserver.com

- 1. Choose Start/Run from your Windows Start menu, type CMD in the Open field and press OK to open a command prompt.
- 2. Type the following command then press [enter] to open a connection between the DawnMailer PC and your SMTP server using port 25:

telnet smtp.myserver.com 25

Telnet connects to your SMTP server and reports its success with a message similar to the following:

220 smtp.myserver.com ESMTP Mail Service Version: [version-number] ready at [today's date & time]

3. Next declare your location by issuing the following:

HELO local.domain.name [enter]

...where *local.domain.name* is supposed to be your exact fully qualified domain name. This is not validated however, so if you are unsure you can type something like

yourfirstname.yourlastname.org

The SMTP server replies by acknowledging your ip address:

250 mail.domain.ext Hello local.domain.name [your ip address], pleased to meet you

4. Declare the email address you wish to send the email from (ensure this is the address you are using as the Sender address in Dawn so you can see whether or not it is accepted), by issuing the following command

MAIL FROM: dawnsender@dawnsenderaddress.org [enter]

If the SMTP Server accepts the sender address, it returns the following message:

250 2.1.0 dawnsender@dawnsenderaddress.org ... Sender ok

5. Declare the recipients email address by typing the following command (it's a good idea to use your own email address here so you can check whether the message is sent):

RCPT TO: myname@myaddress.com [enter]

If the SMTP accepts the recipient address it returns the following message

250 2.1.0 myname@myaddress.com... Recipient ok

- 6. Issue the command, DATA [enter] to start composing the message
- 7. If you want a subject for your email type Subject:-type subject here-then press enter twice.
- 8. Type the body of your message e.g "Test message sent via Telnet"

When you have finished enter a single "." on a line on it's own and press [enter]. The mail server replies with:

250 2.0.0 Message accepted for delivery

9. Issue the command, QUIT [enter] to close the connection.

If you encounter an error at any point you now have a better idea of what the problem is e.g.

501 noreply@medical.org	You have given an unknown email address (or at least an unknown domain)
Sender domain must	for the sender's email
exist	
550 dawn@dawnaddress.com	The mail server has refused to relay mail for you. This may be because:
Relaying Denied	SMTP server does not permit relaying except for certain IP addresses
550 dawn@unknown.com	(this is a common defense against SPAMMING). In this case you may have to talk to the administrator about adding the DawnMailer machine
Not Permitted	to the permitted list.
	 Not using an email address provided by the owner of the server. (This may be the case if you are using a commercial 3rd party provider)



Appendix 1: Diagram - Print Station / Dawn Mailer Printer Mappings

F:\MARKETING\Website\New 4S Website\Content\1. Anticoagulation\1.12 Customer Support Resources\DawnMailerTroubleshootingGuide.doc Version 1.00 Page 23 of 23