



DAWN Reports Webpage

- A web link has been emailed out to all lead users.
- The support team have created many custom reports for sites over the last couple of years.
- We want to share these with you.

<http://www.4s-dawn.com/dawnac/Reports.htm>



DAWN Reports Webpage

Makes life easier - can download and install reports in minutes.

No need to wait for 4S team to create new report.

Page gets updated as and when we create new reports that are not already on the list.



DAWN Reports Webpage

Reports split into two sections:

- List reports, to display lists of patients or information
- Statistics reports, to display counts, sums, figures



DAWN CLINICAL FRAMEWORK

Report name and description of what the report provides is displayed here

Hovering your cursor over the link here displays a preview of what the report looks like

Report Name	Description	Report Preview	Word Instruction Document
% of INRs Between a Certain Value for a Date Range	Reports the % of INRs between certain specified values for a date range.	Hover Here	
% INRs in Range by Target Range	Reports the % INRs in range. This report can be filtered by date range and specific target range.	Hover Here	
% Time in Range with Diagnosis Filter	Reports the % time in range broken down by target range, with the option to filter on a particular diagnosis.	Hover Here	
% Time in, above and below Range by Clinic and Diagnosis	Reports time in, above and below range by clinic and diagnosis, optionally for a date range.	Hover Here	% Time in Above and Below Range By Clinic and Diagnosis

Reports

Report

% of INRs between values

Start date

01/08/2011

End date

30/09/2011

INR Greater Than or Equal to

2

INR Less Than or Equal to

3

[Show report](#)
[Download as XML](#)
[Download as Text](#)

Email address
[Send](#)

Total INRs

% INRs Between Values

26


84.62



DAWN CLINICAL FRAMEWORK

Click on the link here to download or open a Word document. This contains instructions on how to install the report.

Report Name	Description	Report Preview	Word Instruction Document
% of INRs Between a Certain Value for a Date Range	Reports the % of INRs between certain specified values for a date range.	Hover Here	% INRs between Values for a Date Range
% INRs in Range by Target Range	Reports the % INRs in range. This report can be filtered by date range and specific target range.		
% Time in Range with Diagnosis Filter	Reports the % time in range broken down by target range, with the option to filter on a particular diagnosis.		
% Time in, above and below Range by Clinic and Diagnosis	Reports time in, above and below range by clinic and diagnosis optionally for a date range.		





DAWN CLINICAL FRAMEWORK

Report : % INRs Between Values for a Date Range

To Set Up the Report

This report returns the % of INRs between two specific values, for a date range.

1. Press F5 to display the system menu

2. Choose Reports/Report from the Management options to display the Report set up form.

3. Press the New button to create a new report. A blank report form is displayed.

New

4. Enter an appropriate name and description for the report for example: '% INRs between Values'
5. Copy and Paste the following (SQL) text into the top box on the Query tab:

